Product information

Model : BVQ-1100

Product: : BVQ-1100 Real-time PCR Analyzer

Manufacturer

Bioguard Corporation

Address : 6F, No.20-1, Wu-kong 5th Rd., Wugu District, New Taipei City 24890, Taiwan

Tel : +886-2-55916688

Version Information

This manual may be updated without further notice.

Version : 1.0

Issue date : 02.2023

Bioguard Corporation

International Symbols Description

Symbol	Description	Symbol	Description
	Keep away from open flame		Keep dry
	This side up		Handle with Care
	Fragile		Keep away from sunlight

▲ Table 1 Symbol Description

Technical Support and Warranty

Maintenance

This product is guaranteed for one year. Please keep the packing box for transportation, storage, or for use when product needs to be shipped for maintenance.

Bioguard is responsible for the safety, reliability and performance of this product under the following conditions:

- ◆ The analyzer is assembled, upgraded, or maintained by an authorized party.
- ◆ Relevant electrical equipment complies with national standard.
- ◆ This product is used in accordance with this instruction manual.

If the product is in the following circumstances, which is not covered by the warranty, the consumer may need to pay the full amount for maintenance costs.

- ◆ The product is damaged due to the following reasons:
 - Improper operation
 - Improper connection with other devices
 - Accidents
- ◆ Modification of products without authorization
- ◆ Unrecognized product serial number

Technical Support

If technical support is needed, please contact your local distributor.

Please identify the following information:

- ◆ Type of malfunction or damage
- ◆ Product model and serial number

Safety Notice

Warning messages

The terms "danger", "warning", "caution", and "attention" used in the manual indicate hazards and remind users to pay attention to the process of operation according to the degree of hazards. A hazard is a source of potential harm or adverse effects to the environment or user.

	Attention	Provide users with useful information and tips to ensure maximum performance of this product.
	Caution	Refers to a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.
	Warning	Refers to a hazardous situation which, if not avoided, could result in death or serious injury.
	Danger	Refers to an immediate hazard which, if not avoided, will result in death or serious injury.

Safety Information

	Warning	<p>Do not let the lens or image processing unit get wet and avoid the splashing reagents onto the device.</p> <p>Make sure that the ambient temperature during operation is kept at 10 ~ 35° C</p> <p>Please treat the amplification products of this instrument as biohazards and dispose of it properly.</p> <p>In order to safely use this product, the hospital must provide a 110V/220V, 50/60Hz mains power output socket with standard protective grounding.</p> <p>Do not plug or unplug the cable until the entire system is fully switched off.</p>
	Caution	<p>Without the permission or written authorization of the manufacturer, the user shall not disassemble, modify or use various accessories that do not meet the manufacturer's regulations.</p> <p>Please switch off and unplug the device before cleaning it.</p> <p>Do not disassemble the instrument. Doing so may cause an electric shock.</p> <p>Please follow the instruction for cleaning the instrument.</p> <p>In order to prevent damaging the instrument, please carefully follow the instructions described in this manual.</p> <p>This instrument should be regularly calibrated and maintained by a qualified technical personnel.</p> <p>If you have any questions, please contact your local distributor.</p>

Table of contents

Chapter 1 Introduction	6
1.1 Intended use.....	6
1.2 Purpose.....	6
1.3 Product features.....	6
1.4 Introduction of Qmini BVQ-1100 real-time PCR analyzer	7
1.4.1 Components.....	7
1.4.2 Technical specifications and features	8
1.4.2.1 Analyzer specifications	8
1.4.2.2 Technical specifications	8
1.5 Unpacking and shipment.....	9
Chapter 2 Installation	9
2.1 Unpacking.....	9
2.1.1 Packaging contents	9
2.1.2 Unpacking guidelines	10
2.2 Environment.....	10
2.2.1 Operating environment.....	10
2.3 External connection	11
2.4 Installation.....	11
Chapter 3 Getting started.....	12
3.1 Inspection before starting up	12
3.2 Starting the instrument	12
Chapter 4 Software.....	13
4.1 Introduction.....	13
4.2 Function.....	13
4.2.1 Main interface.....	13
4.2.2 Previous data.....	14
4.2.3 Management System	14
4.3 Instrument Settings	15

Table of contents

4.3.1	WiFi settings	15
4.3.2	Date and time settings	18
4.4	Test.....	19
4.4.1	Setting the sample name and test items	20
4.4.2	Testing	22
4.4.3	Checking real-time testing results.....	24
4.5	Previous Data.....	24
4.5.1	Checking previous data	24
4.5.2	Amplification curve.....	26
4.5.3	Checking.....	27
Chapter 5	PC Software	27
5.1	Software Installation	27
5.2	Connecting the instrument	29
5.3	Data update.....	29
5.4	Previous data.....	30
5.4.1	Testing Report.....	30
5.4.2	Checking the amplification curve	32
Chapter 6	Cleaning and Maintenance.....	33
6.1	Cleaning.....	33
6.1.1	Instrument surface cleaning	33
6.1.2	Sample well cleaning.....	33
6.2	Maintenance.....	33
Chapter 7	Troubleshooting	34
7.1	Troubleshooting Guides	34
Chapter 8	Technical support	35

Chapter 1 Introduction

This manual provides general information and operating instructions of the real-time PCR instrument. This chapter provides an overview of hardware, software, and technical support information.

Summary of operating, applications and features:

The analyzer is designed for performing real time polymerase chain reaction (RT-PCR). It supports single or multiplex RT-PCR detection of interested targets in the clinical specimen. A color touchscreen and intuitive software makes the analyzer user friendly. The closed reaction system of the analyzer can reduce contamination and increase productivity. The software can collect and analyze data in real time, enables users to monitor data and print testing reports.

1.1 Intended use

The analyzer is intended to detect DNA/RNA of targeting pathogens using fluorescence-based reagents.

1.2 Purpose

Detecting infectious pathogens of pets.


1.3 Product features

The Real-time PCR analyzer provides rapid and accurate detection of multiple genes. Testing report can be monitored and printed.

1.4 Introduction of Qmini BVQ-1100 real-time PCR analyzer

1.4.1 Components

This product is composed of heating and cooling system, temperature monitoring system, fluorescence detection system, cooling system, and power cord, as shown in Fig1.

▲ Fig. 1 Components of Qmini BVQ-1100 real-time PCR

- | | | | |
|-----------------|----------------|---------------|-----------------|
| 1. Lid switch | 2. Lid | 3. Left panel | 4. Screen |
| 5. Label | 6. Vent | 7. RJ45 | 8. USB port 2 |
| 9. Power switch | 10. Power port | 11. Case feet | 12. Sample well |
| 13. Heating lid | | | |

1.4.2 Technical specifications and features

1.4.2.1 Analyzer specifications

Analyzer	Connectivity
Dimensions: 268*236*188 mm	Network connection : TCP/IP, Ethernet
Weight: 4.5 kg	Operating Conditions
Power Supply	Temperature : 10°C -35°C
Voltage: AC 110V/220V	Humidity : 10%-95%RH
Frequency: 50HZ-60HZ	Atmospheric pressure : Standard atmospheric pressure
	Altitude : 2000 meters
	Storage: Away from moisture
	Noise: 50dB

1.4.2.2 Technical specifications

Thermal parameters

- Temperature accuracy $\pm 0.25^{\circ}\text{C}$
- Temperature uniformity $\pm 0.25^{\circ}\text{C}$
- Temperature control accuracy $\pm 0.1^{\circ}\text{C}$
- Rate of Heating and Cooling $6^{\circ}\text{C} / \text{S}$

Optical parameters

- FAM emission $520\text{nm} \pm 10\text{nm}$
- VIC emission $570\text{nm} \pm 10\text{nm}$
- FAM excitation $495\text{nm} \pm 10\text{nm}$
- VIC excitation $535\text{nm} \pm 10\text{nm}$

Interaction system parameters

- Touch screen : 7" Color touch screen
- USB port : Connect USB flash drive and barcode scanner

Prohibit : The USB port of this instrument cannot be used for charging electronic devices

Note: Before connecting the USB flash drive to the instrument, please format the USB flash drive to FAT32, and confirm that there is enough space (1 GB recommended) in the USB.

1.5 Unpacking and shipment

Do not turn the instrument upside down. After taking out the instrument and accessories, keep the packaging materials for use in shipments for repairs. Check the accessories according to the packing list.

In addition, you should receive all the optional parts and accessories you ordered. If there is anything missing, please contact our local distributor. If damage occurs during transportation, keep the packing materials for inspection and contact your local distributor immediately. According to the packing list provided by the manufacturer, please carefully check if the instrument and accessories are correct, and if there is any damage. Please fill in the product warranty card carefully and send it back to our company within one month of receiving the product.

Chapter 2 Installation

2.1 Unpacking


2.1.1 Packaging contents

When you receive our Qmini BVQ-1100 real-time PCR instrument, please check that it includes the following items.

1.	Qmini BVQ-1100 real-time PCR		1
2.	Pipette	100-1000 ul	1
3.	Pipette	2-20 ul	1
4.	Pipette tip	1 ml	1
5.	Pipette tip	200 ul	1
6.	USB flash drive		1
7.	User guide (simple version)		1
8.	User manual		1
9.	Power cord	International standard, 1.5 meter, black	1
10.	Universal adapter		1
11.	Certificate		1
12.	List of product contents		1

2.1.2 Unpacking guidelines

This product is transported with well-designed carton packaging. The periphery and corners of the instrument are supported and fixed by the inner lining protective foam to prevent collision and vibration during transportation, as shown in Figure 2.

▲ Fig. 2 Packing lining

2.2 Environment

This instrument should be used and kept in an indoor environment.

2.2.1 Operating environment

1. Temperature: 10°C - 35°C
2. Relative humidity: 10% - 90%
3. Keep the instrument on a flat, stable, and dry surface.
4. The instrument should be operated at an altitude of below 2000 meters.
5. Do not place the instrument in a humid or dusty place. The room should be well ventilated and free from corrosive gas or strong magnetic field interference.
6. Do not cover or block the vents in the instrument to prevent it from overheating. The space between the vents of the instrument and the surrounding objects should be at least 30 cm.
7. High temperatures can affect the performance of the instrument. This instrument should avoid direct sunlight and strong light sources. Keep away from heaters, stoves, and other heat sources.

2.3 External connection

Power requirement: This instrument uses a power switch. The operating voltage is 110V/220V, and the frequency is 50/60HZ. Power line is a single phase 3-wire system and should be reliably grounded.

Warning: Electricity needs to be grounded to help safeguard the instrument from any broken circuits!

2.4 Installation

1. The operating environment of the instrument should meet the environmental requirements of 2.2.1 Keep the instrument on a flat, stable, and dry surface.
 2. Please use the power cord provided with the instrument. When connecting, set the power switch to the OFF position.
-

- If the power cord is broken, it must be replaced.
 - Replacement must be replaced with a power cord of the same type and size.
-

Chapter 3 Getting started

3.1 Inspection before starting up

Before turning on the power, check the following :

- 1) The power supply meets the system requirements.
- 2) Make sure the power cord is correctly inserted into the power outlet
- 3) The operating environment of the instrument meets the requirements.


3.2 Starting the instrument

Connect the power cord to the instrument, plug it into a power outlet, and set the power switch to “-”. The following display will appear (Fig 4), and the device will do self-testing.

▲ Fig. 4 Bioguard’s Logo on the screen

Then the main interface appears, as shown in Fig 5.

▲ Fig. 5 Main interface

Chapter 4 Software

4.1 Introduction

The software can help users not only set up and run the tests, but also collect and analyze data. In addition, users can check real-time amplification curve during the PCR amplification.


4.2 Function

Functions of software include:

1. Key in sample names
2. Select testing items
3. Run the test
4. Check the results
5. Check previous data
6. Set up system

4.2.1 Main interface

The main interface will display top bar, middle bar and bottom bar. The upper left corner of the top bar displays the name of the operation interface, and the upper right corner displays the time and date. The middle bar displays the specific operation functions; the bottom bar displays the power button on the left and the LOGO on the right, as shown in Figure 6.

▲ Fig. 6 Main interface

4.2.2 Previous data


Click the “Review”, and enter the previous data interface, as shown in Fig 7. Click “Return” to return to the main interface.

Sample ID	Detection Item	Ct Value	Results	Time
hhh	Feline parvovirus	0	Negative	2022-04-17 16:01:26
hhh	Alphacoronavirus 1	0	Negative	2022-04-17 16:01:26
hhh	Giardia	0	Negative	2022-04-17 16:01:26
hhh	Tritrichomonas foetus	0	Negative	2022-04-17 16:01:26
ggg	H. pylori	0	Negative	2022-04-17 16:01:26
gggg	Leptospira	0	Negative	2022-04-17 16:01:26
gggg	Toxoplasma gondii	0	Negative	2022-04-17 16:01:26
gggg	Babesia canis	0	Negative	2022-04-17 16:01:26

▲ Fig. 7 Previous data interface

4.2.3 Management System

Press the LOGO icon “bioguard” in the lower right corner of main interface for 5 seconds to enter the setting interface. In the management system, you can check device information. Adjust date and time. Select different language. Adjust screen brightness. Import and export data. Set printers. Set IP addresses and enter the maintenance or service interface. Click “Return” to return to the main interface.


▲ Fig. 8 Setting interface

4.3 Instrument Settings


4.3.1 WiFi settings

After entering the setting interface shown in Fig 8, click “Communication Settings”, and enter the network interface shown in Fig 9. Click “Scan WiFi” and enter the interface shown in Fig 10. Click the drop down arrow and select the WiFi as shown in Fig 11. Enter the WiFi password (Fig 12), and click “connect”. When the WiFi is connected, the WiFi logo will appear on the right side of the top frame (Fig 13). If the WiFi is not connected, please check if the password is correct, or place the instrument near the WiFi, and then try to connect again.


After connecting to WiFi, you can click “WiFi Details” to get the WiFi IP address (Fig 14). Please check Chapter 5 for more detailed information regarding the software.


▲ Fig. 9 Network interface


▲ Fig. 10 Scan WiFi


▲ Fig. 11 Searching WiFi

▲ Fig. 12 Enter WiFi password


▲ Fig. 13 Connect WiFi


▲ Fig. 14 Getting WiFi IP

4.3.2 Date and time settings


Click time icon as shown in Fig 8, and enter the date and time setting interface shown in Fig 15. Enter the current date. There are 8 different date formats for selection as shown in Fig 16. You can also enter the current time. When the settings is completed, click “Return”, and a message box will pop up for confirmation of new settings, as shown in Fig 17. When confirming modification, please click “Yes”. If you do not want to keep the settings, click “No”. The date and time are modified from the factory time 2022-11-10 13:46:56 to the current time 2023-02-10 12:00:00 as shown in Fig 18.


▲ Fig. 15 Date and time settings

▲ Fig. 16 Select date formats

▲ Fig. 17 Saving date and time


▲ Fig. 18 Modified date and time

4.4 Test

This section introduces the Qmini BVQ-1100 real-time PCR software and how to complete a PCR test, including sample information settings, data collection, and data management.


4.4.1 Setting the sample name and test items

Turn on the machine and enter the main interface shown in Fig 19, press the "OPEN" button of the instrument to open the lid, put the samples into the sample wells shown in Fig 1, and then click the sample to select the corresponding well in the main interface of the instrument (For example, if sample 1 is placed in the No. 1 well of the instrument, click the No. 1 position button on the main interface), and you will enter the sample editing interface shown in Fig 20.


▲ Fig. 19 Main interface

First, select the species of pet (cat, dog, or exotics). After entering the corresponding interface, click to select the corresponding test item, such as coronavirus (Fig 20).


▲ Fig. 20 Sample editing interface

Then enter the sample ID. For example, we placed the sample of "mxr" in sample well 1, and correspondingly, we enter "mxr" into the sample ID input box (Fig 21)


▲ Fig. 21 Entering sample ID

After entering the sample ID and testing items, click the "Save" button to return to the main interface. At this time, the main interface will display the information entered in the previous step.

▲ Fig 22 Main interface after entering information of sample

Repeat the above steps, and enter other sample information. Up to 8 sample information (8 sample wells) can be entered at the same time, as shown in Fig 23.


▲ Fig. 23 Main interface after entering 8 sample information

4.4.2 Testing

After setting the sample information, click “Run” to start the test, as shown in Fig 24.

The middle part of the functional area displays the amplification curve of real-time PCR, the position of the sample well, and the operating temperature of the instrument respectively. The lower part of the functional area displays the sample well position button, the result viewing button, the countdown timer, and the stop button.


▲ Fig. 24 Instrument running status 1

Click the thermometer area to display the current testing information, as shown in Figure 24, this area displays the amplification curve with indicated color, testing items, and sample well location information. Click the sample well button in the lower part of the selection function area to display the corresponding real time information of the selected sample well.


If the sample information is collected for the sample well, the button of the sample well is selected by default, and the display is blue-green. After the selection is canceled, the color of the button is green, as shown in Figure 25.

If no sample information is collected, the color of the button is white, and the button cannot be selected at this time.

▲ Fig. 25 Instrument running status 2


If you want to terminate the testing before it is completed, please click "Stop". Then, a notice box, "Do you want to stop the test?", will pop up on the interface. Select "Yes" to stop the test, or select "No" to continue the test, as shown in Fig 26.

▲ Fig. 26 Stop the test

4.4.3 Checking real-time testing results

After the testing is completed, it will display the data processed by the software. Click “View Results” to check the results of selected sample well, as shown in Fig 27.


▲ Fig. 27 Results of real time PCR ◦

When the testing is completed, click “Return”. Then, it will go back to main interface.

4.5 Previous Data


4.5.1 Checking previous data

In the main interface, if no sample information is entered, “Review” will be selected, as shown in Fig 28.


▲ Fig. 28 Main interface

Click "Review" to enter the previous data interface. At this time, the function bar displays previous data, page up and down buttons, and the bottom bar displays the "Query" and "Curve" buttons, as well as the current page number and the total page number. The curve button cannot be selected if no sample is selected, as shown in Fig 29.

▲ Fig 29 Previous data interface


After selecting the item for checking, the selected sample is marked with blue-green color. The "curve" button can be selected at this time, as shown in Fig 30.


▲ Fig. 30 Select the item for checking

4.5.3 Checking

In the previous data interface, click the query to enter query interface. Data can be searched based on different query. As shown in Figs 33 and 34, search the required data by entering the measuring date or sample ID

▲ Fig. 33 Query interface

▲ Fig. 34 Results of query


Chapter 5 PC Software

5.1 Software Installation

Find the USB flash drive in the accessory of the instrument. Insert the USB flash drive into the computer with Windows system, open the USB flash drive (Fig 35), copy the compressed package of “Qmini reporting system (use


directly after decompression).7z” to the computer, and decompress it.

Open the decompressed folder, click the application "Qmini reporting system" to enter the PC software interface.

▲ Fig. 35 PC software folder


The PC software is divided into top bar, middle bar, and bottom bar. The top bar shows the software name and version number. The middle part displays software navigation and specific information, and the bottom part displays the instrument connection status. The left navigation in the middle part has system settings and data management, which can set up the connection between PC and instrument, update instrument data, display the operating status of the instrument in real time, and view previous data, etc. When the PCR instrument is not connected, the indicator icon in the bottom bar is red (Fig 36).

▲ Fig. 36 PC software interface

5.2 Connecting the instrument


Make sure that the PCR instrument is turned on, and then obtain the WiFi IP address of the instrument according to step 4.3.1. Click "Set & Connect" on the PC side to enter the interface shown in Fig 37, and enter the WiFi IP address of the instrument in the IP address box (the input in Fig 37 is 192.168.3.15). Click Connect, and then the indicator icon on the bottom turns blue-green. It means that the instrument is connected, as shown in Fig 37. If the connection is abnormal, please close the software. Check if the WiFi environment communicates normally, check if the instrument is connected well, and repeat the operation steps according to the steps in 5.1.

▲ Fig. 37 Instrument connection

5.3 Data update

Update the testing data of the instrument to the PC: click "Data management" to enter the previous data interface (Fig 38). Click the "Update data" button, and the progress bar will display the progress of updating. When "100%" is displayed, it means that the data update is completed. If the update is slow, please move the instrument to an environment with a stronger WiFi signal.

▲ Fig. 38 Data update

Notice :

Stronger signal strength is correlated with higher data transfer speeds.

5.4 Previous data

Click "Data management" to enter previous data interface. Select the item for checking, as shown in Fig 39.

No.	Sample ID	Detection Item	Ct Value	Results	Time
1	5	Chlamydia	29.58	Positive	2022-09-17 12:17:40
1	5	Mycoplasma	28.56	Positive	2022-09-17 12:17:40
2	5	Chlamydia	29.38	Positive	2022-09-17 12:17:40
2	5	Mycoplasma	29.35	Positive	2022-09-17 12:17:40
3	5	Chlamydia	28.93	Positive	2022-09-17 12:17:40
3	5	Mycoplasma	29.11	Positive	2022-09-17 12:17:40
4	5	Chlamydia	29.72	Positive	2022-09-17 12:17:40
4	5	Mycoplasma	29.86	Positive	2022-09-17 12:17:40
5	5	Chlamydia	28.06	Positive	2022-09-17 12:17:40
5	5	Mycoplasma	28.47	Positive	2022-09-17 12:17:40
6	5	Chlamydia	29.45	Positive	2022-09-17 12:17:40
6	5	Mycoplasma	29.25	Positive	2022-09-17 12:17:40
1	5	Chlamydia	29.39	Positive	2022-09-16 10:30:29
1	5	Mycoplasma	29.03	Positive	2022-09-16 10:30:29
2	5	Chlamydia	30.34	Positive	2022-09-16 10:30:29
2	5	Mycoplasma	29.94	Positive	2022-09-16 10:30:29

▲ Fig. 39 Previous data

5.4.1 Test Report


Test report provides:

1. Export PDF reports
2. Print reports

Select the item and click "Report" to enter the report interface (Fig 40). The information editing section is on the left, and the report preview section is on the right. After editing the information, click "Print" to enter the preview interface shown in Fig 41. Then, click the "🖨️" button to complete the printing.

Test Item	Ct Value	Result	Reference
Mycoplasma	28.47	Positive	0-5*10 ⁸ cfu/ml, the result was determined to be positive. Ct>38 or Ct value < 0, the result is judged to be negative.

▲ Fig. 40 Report interface


▲ Fig 41 Report preview

Notice: Make sure the computer is connected to the printer.


To export the report, select the folder to be stored in the storage path section, click "Export Report", and the PDF file of testing report will be generated in the corresponding location (Fig 42).


▲ Fig. 42 Export PDF reports

5.4.2 Checking the amplification curve

Click “Curve” in the main interface to enter amplification curve interface. The amplification curve and sample detection information of selected sample well can be displayed via selecting the number of sample well “1 2 3 4 5 6 7 8” “FAM” and “VIC” checkboxes can distinguish results of dual detection shown in Figs 43 and 44.

▲ Fig. 43 Amplification curve 1

▲ Fig. 44 Amplification curve 2

Chapter 6 Cleaning and Maintenance

6.1 Cleaning

6.1.1 Instrument surface cleaning

The surface of the instrument should be cleaned regularly using a slightly moist soft cloth with water, and the instrument should be wiped dry after cleaning.

6.1.2 Sample well cleaning

1. If the sample well is contaminated with dust or impurities, it will affect PCR amplification and fluorescence detection. Therefore, it should be cleaned regularly, usually once every three months, and can be gently blown with a blower.
2. In order to prevent dust from entering the sample well, the lid should be closed when the instrument is not in use.
3. If any reagent enters the sample well, wipe it clean using a dust-free soft cloth with 70% alcohol.

Notice :

Switch off and unplug the instrument before cleaning.

Do not pour liquids into the reaction module or inside the instrument. Do not scrub the instrument with strong corrosive or organic solvents.

6.2 Maintenance

1. Do not switch the instrument on and off frequently.
2. Do not turn off the power immediately after testing, keep it in standby mode for 10 minutes (the internal fan of the instrument is still working continuously during this period). Turn off the power after the module temperature drops to room temperature.
3. Please use the power cord provided by the manufacturer.
4. Do not use hot water bath or low temperature insulation (such as 4° C) on the instrument.
5. Do not disassemble the instrument. If any problem occurs please contact your local distributor.

Chapter 7 Troubleshooting

7.1 Troubleshooting Guides

No.	Problem	Possible causes	Solution
1.	After booting, the status LED is normal, but the screen is black	Abnormal connection	Check if the circuit is stable, if the voltage is normal, and if the power interface is tightly plugged in
2.	Touch screen is unresponsive.	Faulty screen or poor cable contact	Contact maintenance personnel.
3.	Fan does not work	Internal cable problem or fan failure	Shut down the instrument, and contact maintenance personnel.
4.	There is no response when the USB flash drive is inserted into the USB port.	The format of USB flash drive is not compatible.	Check if the format of USB flash drive is FAT32.
5.	Wi-Fi network is not detected	The wireless frequency of the WiFi network does not match	Check if the network is 2.4G WiFi network
6.	Network error, please check the network	The device is not connected to the Internet, or the connected network signal is unstable, or the router is faulty	Check the network environment
7.	When uploading data, it prompts "unknown error".	The device is not connected to the Internet, or the connected network signal is unstable, or the router is faulty	Check the network environment
8.	Heating abnormal fault prompt	Abnormal heating	Contact maintenance personnel
9.	Fault prompt for abnormal temperature drop of the instrument	Abnormal cooling of the instrument	Contact maintenance personnel
10.	Heating lid sensor abnormal fault prompt	Abnormal heating lid sensor	Contact maintenance personnel
11.	Abnormal fault prompt of the heat sink sensor	Abnormal heat sink sensor	Contact maintenance personnel
12.	Heat sink over temperature fault prompt	Abnormal heat sink	Contact maintenance personnel

Chapter 8 Technical support

Bioguard guarantees that the Qmini BVQ-1100 real-time PCR analyzer has passed a comprehensive test and meets the requirements in the instruction manual.

The use of this instrument must comply with the instructions and safety warnings given in this user manual, otherwise it will not be covered by the guarantee.

Software: Any software related to this product is provided to customers free of charge as a service, and the software is a necessary tool to run the instrument.

Responsibility: Bioguard is not responsible for direct, indirect, or incidental damage caused by failure to comply with the instructions given in the user manual or incorrect use of the Qmini BVQ-1100 real-time PCR instrument. Only our technicians or their agencies can inspect or provide any parts of the instrument. Bioguard is not responsible for any direct, indirect or incidental damage caused by the user's random disassembly and replacement of parts.

The testing results are for in vitro veterinary diagnosis only, and cannot be used as the full or direct basis for decision-making on clinical medication, surgery, and prognosis. The veterinarian should combine the judgment of the clinical symptoms of the pet with other diagnostic results.

Factors may affect the results during PCR testing include:

Environmental factors such as aerosol pollution of the previous testing items, unstable voltage, etc. Human operation factors such as extraction errors, failure to strictly follow the requirements of the instruction manual, etc.

Equipment factors such as vibration, improper use, etc.

Kit storage factors such as high temperature, high humidity, strong direct light, unqualified reagent quality, etc.

Test results may appear false negative, or false positive, or inconsistent with clinical symptoms due to the factors mentioned above. The test results are for in vitro veterinary diagnosis only. Please consult all other clinical symptoms and use other diagnostic methods, to ensure an accurate diagnosis.

Before using the Qmini BVQ-1100 real-time PCR analyzer, please read this manual carefully to ensure that you can get the best testing results.


If you have any questions regarding the testing process, please contact our technical support.

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The screenshot shows the 'History Data' interface. At the top, there is a header bar with a '記錄' (Record) button, signal strength and Wi-Fi icons, and a back arrow with the timestamp '2022-11-10 13:46:56'. Below the header is a table with the following data:

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Below the table, there are navigation buttons: '查詢' (Search), '匯出報告' (Export Report), and '擴增曲線' (Expand Curve). The page number '1/6' is shown at the bottom left, and the 'bioguard' logo is at the bottom right.

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