# Introduction

Thank you for purchasing Qmini BVG-1100 Nucleic Acid Extractor. Before using the instrument, please read this user guide carefully in order to use the product correctly. Please place this manual nearby the instrument for a timely reading if necessary

Product: Qmini BVG-1100 Nucleic Acid Extractor

Model: BVG-1100

## Version Information

This manual may be updated without further notice.

Version: 1.0

Issue date: 01.2023

## Statement

Bioguard does not make any form of guarantee for this information, including (but not limited to) the implied guarantee of merchantability or fitness for a specific purpose.

Bioguard is only responsible for the safety, reliability and performance of the instrument under the following circumstances:

- ◆ The instrument is assembled, upgraded, or maintained by authorized personnel.
- ◆ Relevant electrical equipment complies with national standards.
- ◆ Product is used in accordance with this instruction manual.

Bioguard cannot and does not assume responsibility for the safety, reliability, and performance of the instrument under the following circumstances:

- ◆ Parts are disassembled and readjusted.
- ◆ The user manual was not followed for correct operation.

# Warranty and Maintenance

The warranty period of the product is 1 year and 1 year warranty for major accessories. The warranty period starts from the "installation date" filled in on the "Equipment Warranty Card" attached to the product, and the "Equipment Warranty Card" is the only certificate for calculating the warranty period. To protect your rights, please ask the installer to return the second copy of the "Equipment Warranty Card" to our company within 30 days from the date of installation. If the "Equipment Warranty Card" corresponding to the product you purchased fails to return to the company on time, the warranty period will start on the 45th day from the "delivery date" marked on the equipment box.

During the warranty period, the product can enjoy free after-sales service. However, if the product is in the following circumstances, which is not covered by the warranty, the consumer may need to pay the full amount of maintenance costs.

- ◆ Damage caused by human factors
- ◆ Improper operation
- ◆ The grid voltage exceeds the specified range of this instrument.
- ◆ Natural disasters
- ◆ Replacement or use of parts, accessories, and consumables not approved by Bioguard, or maintenance by a party not authorized by Bioguard
- ◆ Other failures not caused by the product itself
- ◆ After the warranty period expires, Bioguard can continue to provide maintenance services at a charge. If you do not pay or delay the payment of the maintenance service fee, the company will temporarily suspend the maintenance service until payment is made.

# Technical Support

If technical support is needed, please contact your local distributor, or contact the customer service department of Bioguard.

Please provide the following information for the technical support:

- ◆ Type of malfunction
- ◆ Product model and serial number


# Safety Notice

## Warning messages

The terms "danger", "warning", "caution", and "attention" used in the manual indicate hazards and remind users to pay attention to the process of operation according to the degree of hazards. A hazard is a source of potential harm or adverse effects to the environment or user.

	Attention	Provide users with useful information and tips to ensure maximum performance of this product.
	Caution	Refers to a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.
	Warning	Refers to a hazardous situation which, if not avoided, could result in death or serious injury.
	Danger	Refers to an immediate hazard which, if not avoided, will result in death or serious injury.

## Safety Information

	Caution	The instrument must be used with special consumables.
	Caution	Do not interfere the operation of the instrument or move the stationary parts of the instrument during operation.
	Warning	Use a reliable, grounded power outlet.
	Warning	Without the manufacturer's permission or written authorization, the user should not disassemble, modify, or use various accessories that do not comply with the manufacturer's regulations.
	Warning	Before cleaning or disinfecting the instrument, the instrument must be turned off and the power cord must be unplugged.
	Warning	Do not disassemble the instrument to prevent electric shock.
	Attention	Please read the instructions carefully before use.

# Table of contents

Introduction .....	1
Version Information .....	1
Statement.....	1
Warranty and Maintenance .....	2
Technical Support.....	2
Safety Notice .....	3
Chapter 1 System Introduction .....	5
1.1 Introduction.....	5
1.2 Product features.....	5
Chapter 2 Instrument Introduction.....	6
2.1 Components .....	6
2.2 Specifications .....	7
2.3 Operation screen.....	7
2.3.1 Indicator lights .....	7
2.3.2 Key function description.....	8
Chapter 3 Getting Started.....	8
3.1 Packaging contents.....	8
3.2 Operating environment .....	8
3.3 Operation.....	9
Chapter 4 Cleaning and Maintenance .....	10
Chapter 5 Troubleshooting.....	10

# Chapter 1 System Introduction

This manual provides general information and operating instructions for the nucleic acid extractor.

## 1.1 Introduction

Qmini BVG-1100 Nucleic Acid Extractor is an automated nucleic acid extraction equipment used for sample pretreatment for PCR detection.

## 1.2 Product features

**Fast extraction:** 2 samples can be processed at the same time, in as little as 8 minutes.

**Safe and reliable:** The instrument runs automatically without manual operations. Reagent and magnetic cover are clearly marked for easy identification to ensure proper usage and to avoid leakage, ensuring high reliability in the test.


**Contamination control:** Small and closed extraction chamber design with built-in ultraviolet light can effectively remove potential nucleic acid contamination.

**User-friendly:** Built-in operation program with one-button start.


# Chapter 2 Instrument Introduction

## 2.1 Components

Qmini BVG-1100 is a 2-channel fully automatic nucleic acid extractor, which needs to be equipped with special consumables (reagent kit and magnetic cover) for operation.

▲ Fig 1 Appearance of Qmini BVG-1100 Nucleic Acid Extractor

▲ Fig. 2 Qmini BVG-1100 Diagram of dedicated consumables

The magnetic cover is placed in the position as shown on the kit. The instrument can automatically grasp it before the extraction starts and automatically return it to the original position after the extraction is completed. The reagents are preset in the kit and the sample is added to the lysis well before extraction.

## 2.2 Specifications

Methodology: Magnetic bead

### Specifications


Dimensions: 226\*207\*205 mm  
 Weight: 3.7 Kg  
 Power supply: AC 110V/220V 50/60HZ  
 Storage: Keep away from moisture  
 Difference between channels :  $CV \leq 5\%$

### Operation environment

Temperature: 10°C - 40°C  
 Relative moisture: 10% - 90%RH  
 Altitude: 2000 meters  
 Noise: <70dB  
 Magnetic bead recovery : >92%  
 Extraction time : <8min  
 Oscillation : 3HZ  
 Contamination control : UV sterilization

## 2.3 Operation screen

Please refer to the figure below for the key functions and indicators on the operation screen.

▲ Fig. 3 Operation screen

### 2.3.1 Indicator lights

#### Consumables detection indication:

The indicator light has three states: off, red, and green.

Off: Unloaded

Red: Kit detected

Green : Magnetic cover detected

#### Operating status indication :

Standby: Yellow light flashes

Running: Green light on

Disinfection: Blue light on

## 2.3.2 Key function description

**In/out:** Control the reagent kit carrier between the out-of-cabin and in-cabin state. Controlled exit is only valid in the standby state of the instrument and is invalid in the running or disinfection state.

**Running:** Control the operation of the extractor. Press to start the extraction program. It is only valid when the reagent kit is in the cabin and the instrument is in the standby state. It is invalid during the running or disinfection state of the instrument. Long press to forcibly stop the operation.

**Disinfection:** Turn the UV sterilization on or off. Press to turn on the UV. It is only valid when the reagent kit is in the cabin and the instrument is in the standby state. It is invalid in the running state. Press it again to turn off the UV manually, or it will automatically turn off after 15 minutes.

# Chapter 3 Getting Started

## 3.1 Packaging contents

When you receive our Qmini BVG-1100 Nucleic acid extractor, please check that it includes the following items. If you find any missing or damaged items, please contact the distributor or sales representative.

1.	Qmini BVG-1100 Nucleic Acid Extractor	1
2.	Power cord	1
3.	Universal adapter	1
4.	User manual	1
5.	Product warranty	1
6.	Certificate	1
7.	Extraction kit holder	1
8.	PCR rack	1
9.	List of packaging contents	1

## 3.2 Operating environment

1. Temperature: 10°C - 40°C
2. Relative humidity: 10% - 90%
3. This instrument should be used and kept in an indoor environment. Keep the instrument on a flat, and stable surface.
4. Do not place the instrument in a humid or dusty area. The room should be well ventilated and free from corrosive gas or strong magnetic field interference.

## 3.3 Operation

For the first time user, please follow the instructions. Connect the power supply according to the instructions and turn on the power switch. The instrument will immediately reset to its reference point and enter the standby state.

### 1. Prepare Consumables

Take the new reagent kit, shake it well, and place it in the reagent kit carrier. Then, tear off the sealing film, add 200 $\mu$ L of sample to the **LB** well, and place the magnetic cover in its well as shown in Fig 2.

### 2. Place Consumables

Press the key and the carrier will extend out. Then, place the reagent kit into the carrier and confirm that it is properly in place (the 4 supports on the side of the kit will be aligned in the carrier), and note that the magnetic cover should be in the indicated position.

After the kit is placed, the consumable indicator light on the operation screen will change from off to red.

### 3. Extraction :

Press the enter/exit key and wait for the carrier to retreat into the cabin. Then, press the run key.

The extraction program will start to run and the panel status indicator will turn from standby to running. Within 5 seconds, the consumables indicator turns from red to green, indicating that the magnetic cover has been successfully picked up and that the extraction process has begun. After about 8 minutes, the buzzer will sound, indicating that the extraction was completed. The panel status indicator will change from running to standby.

### 4. Eject consumables

Press the enter/exit key and wait for the carrier to extend out. Then, take out the reagent kit and place it on the external carrier. Press the enter/exit key again to retract the carrier. Use a pipette to draw the liquid in the **EB** well for PCR testing.

### 5. Disinfection

Press the disinfection button while in the standby state. The UV lamp in the instrument will be turned on and it will automatically turn off after 15 minutes.

### 6. Turning off the instrument

After sterilization, directly turn off the power of the instrument with the switch.

## Chapter 4 Cleaning and Maintenance

The instrument is a closed structure. If there is no special situation, wiping with alcohol will be sufficient.

Surface cleaning: Wipe with a clean, dust-free cloth and 75% alcohol.

Internal cleaning: Turn on the ultraviolet lamp for disinfection or wipe it with 75% alcohol after extending the reagent kit carrier.

	Danger:	Switch off and unplug the instrument before cleaning.
	Caution:	Do not disassemble the instrument without permission, it may cause damage to the instrument
	Caution:	If the instrument overturns during the operation, causing the reagents to flow out or the consumables to drop, please contact technical support.

## Chapter 5 Troubleshooting

No.	Error message	Possible causes	Solution
1.	After the instrument starts running, the buzzer sounds and the instrument stops running.	The magnetic cover is not in the kit.	Place magnetic cover in the kit as required.
2.	The reagent kit shakes severely when the instrument is running	The magnetic cover is deformed and bent or the kit is not placed in the proper position.	Stop the process and check the magnetic cover and kit position.
3.	No response after pressing the Exit or Run button.	Instrument may be in the disinfection state.	Check and turn off disinfection mode manually.
4.	The instrument cannot operate, and the green light of the consumable detection indicator is on. The buzzer sounds after pressing the run button.	There is a magnetic cover left on the magnetic cover holder.	Load an empty reagent box, then press and hold the run button to perform the demagnetization action. Then, remove the leftover magnetic cover.

If there are other abnormalities, please contact technical support.

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