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

miniCHEM Veterinary Coagulation Test Disc

For Veterinary Use Only

miniCHEM veterinary coagulation test disc is used in conjunction with miniCHEM Veterinary Coagulation and Chemistry Combo Analyzer to quantitatively detect various coagulation function of dogs and cats. It is of great significance for the auxiliary diagnosis of related diseases.

【Packaging】

5 disc/Box
10 disc/Box


【Specimen】

Plasma or whole blood anticoagulated with sodium citrate (blue top tube)


NOTE: Do not use EDTA, heparin or oxalates as anticoagulants. When using a blue top tube, enough blood should be collected up to label mark of the tube.

【Specimen collection】

- The required sample volume for each test is 180 μL ; when using serum or plasma as the sample, please draw 2-3 times the amount of whole blood to ensure a sufficient sample volume.
- Sample collection, storage and testing must use plastic products or siliconized glass products.
- Samples should avoid hemolysis or contamination caused by tissue fluid.
- When using sodium citrate blood collection tubes, enough blood should be collected to the mark on the blood collection tube to maintain the correct ratio of sample to anticoagulant.
- When mixing blood with anticoagulant, invert gently at least five times to mix thoroughly, do not shake vigorously.
- Please use the samples within 60 minutes of sample collection. If it cannot be tested in time, the plasma should be separated and stored at 2 to 8°C, and should not exceed 6 hours.

Take refrigerated disc out and warm it for 10 minutes at room temperature. Then, take the disc out of the sealed bag.


Fix the reagent tray on the tray rack
* Align the groove of the reagent disc with the groove of the tray rack and gently press to fix the reagent disc

- Plasma: Freshly collected venous blood is anticoagulated with 0.109mol/L Sodium Citrate 9 : 1, mix thoroughly, and centrifuge at 3000 rpm for 15 minutes to collect plasma.
- Whole blood : Freshly collected venous blood is anticoagulated with 0.109mol/L Sodium Citrate 9 : 1.

【Testing steps】

1. The reagent disc should be stored at 2~ 8°C. Before use, please take it out and warm it to room temperature for 10 minutes. * The reagent disc is individually sealed and must be used within 30 minutes after opening.
2. Take the reagent disc out of the aluminum bag, and fix the reagent disc on the tray rack. Place it on a horizontal table, use a 200 μL pipette and set the volume at 180 μL . Take 180 μL sample and inject the sample into the reagent disc through the "sample" injection hole (red circle). Take 180 μL diluent and inject the diluent through the "Diluent" injection hole (white circle) into the disc. * Do not touch the tip of the dropper when adding the sample, and avoid adding air bubbles.
3. Open the lid of the instrument, fix the reagent disc horizontally on the tray rack, close the lid and click the start button in the middle to start the test.

Warning: Please wear dust-free gloves for operation. When placing the test disc, please hold the edge of the disc and do not touch other parts of the disc.

【Storage and shelf life】

- Storage temperature: 2°C-8°C; shelf life: 12 months after production
- Do not store in places exposed to direct sunlight or in environments exceeding 30 °C.

【Precautions】

- This product is for in vitro veterinary diagnosis use only.
- Please pay attention to the shelf life before use, and do not use the test disc beyond the shelf life. The test disc is a disposable consumable, please do not reuse it.
- The used test disc, anticoagulation collection tube and dropper tip are all disposable products. Please do not reuse them. After use, please discard them according to your country's biomedical waste disposal regulations.