Pet Pathogen Nucleic Acid Extraction Kit

USER MANUAL

【Composition and specifications】

Nucleic acid extraction reagent 16 units

magnetic cover 16 units

User manual 1 piece

【Storage and shelf life】

Storage temperature 2°C - 8 °C; shelf life: 24 months, please use it within the expiration date.

【Operating environment】

Recommended operating temperature is 20-25 °C.

【Specimens】

Pleural effusion, ascites, stool swabs, anal swabs, eye/nasopharyngeal swabs, EDTA anticoagulant blood, urine.

【Sampling requirements】

Liquid sample: (including blood sample, abdominal effusion, pleural effusion or other liquid samples, etc.)

Take 200 µL of liquid sample into preservation buffer, and mix well. Then, transfer 200 µL of the mixture into the LB well for lysis. For urine sample, directly take 200 µL and add to LB well for lysis.

Swab sample:

1. Fecal swab: take proper amount with a swab
2. Anal swab: moisten the swab with dilution buffer, and then wipe specimen
3. Eye/nasopharyngeal swab: swab under the eyelid and fully wipe to collect eye swab samples; swab oral and nasal discharge properly to collect nasopharyngeal swab samples
4. After swab sample collection, quickly snap the swab handle and place in the preservation buffer, and then shake it to fully dissolve the pathogen on the swab into the preservation buffer.

It is recommended that nucleic acid extraction and detection should be performed immediately after sample collection. If storage needed, samples can be stored at 4°C temporally. Samples need to be refrigerated during transportation.


【Preparation before testing】

1. Please read the instruction manual carefully before testing, be familiar with each step, and strictly follow the requirements when using the kit.
2. Wear disposable gloves and masks, operate under the condition meeting the requirement of PCR testing environment.

【Operation】

Nucleic acid extraction : (using nucleic acid extractor)

1. Take out the extraction cartridges, as needed. Tear off the top seal according to the instruction. The well positions of the rack as shown in below picture.

2. Use the pipette to take 200 µL of sample and add it to the LB well.
3. Take out a magnetic cover and place it in the corresponding position of the cartridge.
4. Press "In/Out ▼" button of the instrument, put the rack in the holder of the automatic extractor, press "In/Out ▼" button again. The rack will be sent into the instrument and the door will close.

Note: The cartridge needs to be placed in the holder according to the direction displayed on the side label and should lie to the bottom of the holder. If the cartridge is located correctly, the indicator light corresponding to the channel position will show red light, otherwise the light will not be on.

5. Press "operate ►" button to start the nucleic acid extraction. After the instrument finish the operation and play a prompt sound, press "In/Out ▼" button and take out the rack.

Note: Extracted nucleic acid samples are recommended to be tested immediately; If not in immediate use, please seal them with stickers and store at -20°C or -80°C.

After nucleic acid extraction, Please follow the PCR amplification instructions for PCR amplification

【Limitations】

1. Aerosol contamination of amplification products can easily cause false positive results. The testing laboratory should be set up strictly following the requirements of the PCR testing laboratory.
2. A negative result cannot completely exclude the possibility of pathogen infection and should be combined with clinical observations and additional testing for interpretation.

【Product quality indicators】

1. Sensitivity: The product detection limit of the kits is 1000 copies/ml.
2. Specificity: The test kits do not detect the cross-reaction of the pathogen samples.
3. Precision: When a strong positive sample and a weak positive sample are repeated testing for 10 times in a row separately, the CV values of their Ct values show less than 5%.

【Precautions】

1. Product quality inspection: Before using PCR amplification reagent, unpack to check if the lyophilized reagent on the tube bottom is normal (white, clumped). If liquefied, it cannot be used. Otherwise, it will affect the PCR results.
2. Pipette use: When drawing 20µL of nucleic acid supernatant as the template for PCR amplification, only depress the plunger of the pipette to the first stop. Do not press to the bottom, otherwise the sample volume will exceed 20µL which will affect the results.
3. Sample information setting: Make sure that the PCR tubes are placed in the same sample well position as set in the instrument. For example, if the PCR tube is placed in the #1 well of the sample plate, then select the corresponding #1 position in the interface. After setting the sample name and the test item, press the Run icon to start the PCR amplification process.

Bioguard Corporation
4F., No. 25, Wugong 5th Rd., Xinzhuang Dist.,
New Taipei City 248020, Taiwan
Tel: +886-2-55916688
<https://www.bioguardlabs.com>

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